

先进的生产设备 Advanced production equipment

全自动浸胶生产线，产品一致性强，生产效率高。
Fully automatic dipping treatment production line improves product consistency and production efficiency.



优秀的生物相容性 Excellent biological compatibility

产品在去离子水中沥滤处理，有效去除有害化学物质，提高生物相容性。

Products are leached in deionized water to remove harmful chemical and improve product biological compatibility.



严格的品控 Strict quality control system

每一支产品都进行了性能及球囊检测，保证每支产品品质。


Every single product is tested to ensure the product quality



单腔乳胶导尿管：用于临时导尿
1-way latex foley catheter: for temporary catheterization
Balloon volume: w/o balloon
Length: 400mm Size : 6Fr-26Fr



乳胶吸痰管：供临床吸痰用
latex Suction catheter: for Suction
Balloon volume: w/o balloon
Length: 400mm Size : 6Fr-26Fr



单腔球囊型乳胶导尿管（扩宫管）：用于子宫扩宫
1-way balloon latex catheter (Uterus Expansion catheter):
For the Expansion of uterus
Balloon volume: 5/10 ml/cc
Length: 400mm Size : 12Fr-16Fr



双腔儿童型乳胶导尿管：用于儿童留置导尿或止血导尿
 2-way latex foley catheter for children: for children catheterization or hemostasis catheterization
 Balloon volume : 3 /5 ml/cc
Length: 300mm Size : 6Fr-10Fr



双腔成人型乳胶导尿管：用于留置导尿或止血导尿
 2-way adult latex foley catheter: for Catheterization or hemostatic catheterization
 Balloon volume: 5/10/15/30 ml/cc
Length: 400mm Size : 12Fr-26Fr



双腔女用型乳胶导尿管：用于女性留置导尿或止血导尿
 2-way female latex foley catheter: for women Catheterization or hemostatic catheterization
 Balloon volume: 5/10/15/30 ml/cc
Length:270mm Size : 12Fr-26Fr



双腔球囊型乳胶导尿管（子宫造影通水管）：
 用于子宫通水、注药、造影时一次性使用
 Uterine radiography catheter : for Uterine radiography
 For the Expansion of uterus
 Balloon volume:5ml/cc
Length:400mm Size : 12Fr-14Fr



三腔标准型乳胶导尿管：用于留置或止血导尿，膀胱滴注药液
 3-way standard latex foley catheter: for Catheterization or Hemostatic catheterization, Bladder infusion drug solution
 For the Expansion of uterus
 Balloon volume: 5/15/30/50 ml/cc
Length: 400mm Size : 12Fr-26Fr



三腔双囊型乳胶导尿管：用于导尿或止血导尿，膀胱滴注药液，前列腺灌注及手术
 3-way Double balloon latex catheter: for Catheterization or Hemostatic catheterization, Bladder infusion drug solution, Prostate perfusion and surgery
 Balloon volume: 5/15/30/50 ml/cc
Length: 400mm Size : 12Fr-26Fr